

AEU-7000

Implant & Oral Surgery System



Features

- Powerful, brushless 40,000 rpm autoclavable micromotor
- Variable handpiece ratio selection: 1:5 and 1:2 increaser, 1:1, and 20:1 reduction E-type handpieces
- Adjustable torque up to 60 Ncm enables a single, high-efficiency 20:1 reduction handpiece (AHP-85MB) to be used for the complete implant procedure
- Advanced calibration technology for accurate speed and torque performance
- Selectable auto-stop function when desired torque setting is reached
- Contemporary design with easily cleaned surfaces
- Fully integrated, easy-load irrigation pump
- Large, bright, easy-to-read display
- Intuitive keypad with adjustable handpiece ratio, speed, torque, and irrigation flow settings
- Six programmable preset buttons for storing and labeling custom settings in memory
- Autoclavable motor holder may be free-standing or mounted to either side of the console
- Upgradeable software for a longer return on investment
- Includes AE-7PM On/Off Foot Control and irrigation tubing set
- Standard tubing set accepts standard irrigation bags; optional tubing set accommodates rigid irrigation bottles
- Auto-sensing global voltage compatibility
- FDA, UL, and CE compliant
- Made in the USA

Foot Switch Upgrade



Optional AE-70V Variable Speed Foot Control with Pump On/Off, Micromotor Direction, Preset Selection, and Torque Cycle buttons



1-800-426-5913 or 425-487-3157
P.O. Box 1548 • Woodinville, WA 98072
www.aseptico.com • info@aseptico.com